

TeleRPM™

TeleRPM AnyHub 2022

The world's first Hub designed for telehealth



Out-of-the-box




No App and
Additional Set up



4G Cellular



Connect with
Transtek Default
BLE Devices



Product Name	TeleRPM AnyHub 2022
Dimensions	100mm in diameter, 48mm in height
Color	Snow White
Accessory	5V/2A power adapter, USB-TypeC cable
Button/Control	Sync Button, Reset Button
Cellular	2G/4G(LTE-Cat-M1), Global Roaming Network
Bluetooth	Bluetooth 5.0
Bluetooth Range	up to 30 meters
Bluetooth Connection	up to 10 simultaneous connection



Contact us

inquiry@transtekcorp.com
www.transtekcorp.com

LinkedIn



Follow us on : Transtek Medical