

TeleRPM™

# TeleRPM BPM Gen 1

Connected Blood Pressure Monitor



Fits 9" -16.5"  
cuff size



One-click  
Operation



4G Cellular



Clinically  
Validation Accuracy



Product Name	TeleRPM BPM Gen 1
Dimensions	Approx.159.2 x 121.5 x 68.1mm
Cuff Size	22-42cm
LCD Display	V.A. 78 x 92 mm
Gross Weight	Approx.393g (Excluding the batteries and cuff)
Measurement Method	Oscillometric technology
Rated Cuff Pressure	0~299 mmHg(0.0~39.9kPa)
Measurement Range	SYS: 60~230 mmHg(8.0~30.7kPa) DIA: 40~130 mmHg(5.3~17.3kPa) Pulse value: 40-199 beats/min
Accuracy	Pressure: ±3 mmHg; Pulse: ±5%
Normal Working Condition	Temperature: +5℃ to+40℃ Relative humidity: 15% to 90%
Storage & Transportation Condition	Temperature: -20℃ to+60℃ Relative humidity: ≤93%
Power Supply	Battery powered mode: 6V DC 4×AA batteries AC adapter powered mode: 6V 1A
Ingress Protection	IP21



LinkedIn

## Contact us

[inquiry@transtekcorp.com](mailto:inquiry@transtekcorp.com)  
[www.transtekcorp.com](http://www.transtekcorp.com)



Follow us on : Transtek Medical