

# TeleRPM™

## TeleRPM Scale (BLE) 2022

Your best weight management partner



441lb Capacity



3.22" LED  
Display



Step-on  
Technology



Bluetooth 5.0

# TeleRPM Scale(BLE) 2022

TeleRPM™



Product Name	TeleRPM Scale(BLE) 2022
Display Size	LED display 80 x 40mm
Scale Size	300 x 300 x 25.1mm
Device Weight	Approximately 1.67kg
Measurement Units	kg/lb
Measurement Range	5 ~200 kg/11 ~ 440 lb
Division	5~100kg, 0.05kg/0.2lb; 100.1~200kg, 0.1kg/0.2lb
Accuracy	Temperature requirement: 25±5°C 50±0.2kg; 100±0.3kg; 150±0.4kg; 200±0.5kg;
Power Supply	6.0V 4xAAA
Power-ON	Press-On
Power-OFF	Auto-off (About 15 seconds while showing 0.0; About 15 seconds after the weight data is locked)
Operating Conditions	Temperature: 5°C to 40°C; Relative humidity: ≤90% RH; Atmospheric pressure: 86kPa to 106kPa
Storage Conditions	Temperature: -20°C to 60°C; Relative humidity: ≤90% RH; Atmospheric pressure: 50kPa to 106kPa
Ingress Protection	IP21



LinkedIn

## Contact us

[inquiry@transtekcorp.com](mailto:inquiry@transtekcorp.com)  
[www.transtekcorp.com](http://www.transtekcorp.com)



Follow us on : Transtek Medical