

TeleRPM™

TeleRPM Scale Gen 1

Connected Weight Scale



551lb Capacity



Automatically
upload data



Step-on
Technology



4G Cellular



Product Name	TeleRPM Scale Gen 1
Scale Size	360 x 360 x 37.2mm
Display Size	85x53mm
Device Weight	6.7lb
Measurement Units	Pound/Kilogram/Stone
Measurement Range	11 ~ 551lb/5 ~ 250kg/0st:11lb to 39st:51b
Division	0.1kg/0.2lb
Accuracy	50kg±0.2kg, 100kg±0.3kg, 150kg±0.4kg, 200kg±0.5kg, 250kg±0.8kg
Power Supply	6V(4x1.5V AA alkaline batteries)
Power-ON	Step on
Power-OFF	Automatic off (10s after showing 0.0)
Operating Conditions	Temperature: 5 °C to 40 °C
Storage Conditions	Temperature: -20 °C to 60 °C
Ingress Protection	IP22



LinkedIn

Contact us

inquiry@transtekcorp.com
www.transtekcorp.com



Follow us on : Transtek Medical