

TeleRPM™

TeleRPM Scale Gen 2

A cellular-enabled scale



One-step
Operation



Only 4.76lb



4G Cellular



Anti-rollover
Design

TeleRPM BPM Gen 2

TeleRPM™



| Product Name | TeleRPM Scale Gen 2 |
|----------------------|---|
| Display Size | LCD display 85x53mm |
| Scale Size | 320x320x28.5mm with tempered glass platform |
| Device Weight | 2.16kg (Net Weight) |
| Measurement Units | kg/lb/st.lb |
| Measurement Range | 5 ~ 200kg(441lb) |
| Division | 0.1kg/0.2lb |
| Accuracy | 50kg±0.2kg, 100kg±0.3kg, 150kg±0.4kg, 200kg±0.5kg |
| Power Supply | 6V 4xAA Batteries |
| Power-on | Step on |
| Power-off | Automatic off (10s after showing 0.0) |
| Operating Conditions | Temperature: 5°C to 40°C |
| Storage Conditions | Temperature: -20°C to 60°C |
| Ingress Protection | IP21 |



Contact us

inquiry@transtekcorp.com
www.transtekcorp.com

LinkedIn



Follow us on : Transtek Medical